


Low speed handpiece bearing 3.967×7.938×3.175mm, smooth, steel balls



Bearing for low speed handpiece, dimensions 3.967×7.938×3.175mm, groove, steel cage, steel balls, bearing closed on both sides, Tecnomed Italia brand.

Shipping expected in approximately 15 days

SKU: 701TM0194

Categories: Internal diameter 4-5

Description

Low speed handpiece bearing, 3.967×7.938×3.175mm, groove type, Tecnomed Italia, spare part for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: height 3.175mm, internal diameter 3.967mm, external diameter 7.938mm
- Groove bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty