

## Bearing, low speed, 3,175x7x2 steel,Z



Bearing for low speed handpiece, dimensions 3.175x7x2mm, smooth, steel cage, ceramic balls, closed on one side, Tecnomed Italia brand.

### Available

**SKU:** CM333

**Categories:** Internal diameter 2-3

## Description

---

Low speed handpiece bearing, 3.175x7x2mm, smooth type, Tecnomed Italia, spare part for dental handpieces, with steel cage and ceramic balls.

### Technical features

- Dimensions: height 2mm, internal diameter 3.175mm, external diameter 7mm
- Smooth bearing
- Steel cage
- Ceramic balls
- Radial contact
- ABEC 7 class
- Bearing closed on one side
- Tecnomed Italia brand
- 12-month warranty