



Bearing, low speed, 4.763X7.938X3.175,steel,ZZ,LI,Bienair MX compatible

Bearing for low speed handpiece, dimensions 4.763x7.938x3.175mm, groove, steel cage, steel balls, bearing closed on both sides, Tecnomed Italia brand.

Shipping expected in approximately 15 days

SKU: 701TM0193

Categories: Internal diameter 4-5

Description

Low speed handpiece bearing, 4.763x7.938x3.175mm, groove type, Tecnomed Italia, spare part for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: height 3.175mm, internal diameter 4.763mm, external diameter 7.938mm
- Groove bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty