

Bearing, low speed, 4x7x1,6 steel



Bearing for low speed handpiece, dimensions 4x7x1.6mm, smooth, steel cage, steel balls, Myonic brand, compatible with Sirona and W&H.

Available

SKU: CN56W

Categories: Internal diameter 4-5

Description

Low speed handpiece bearing, 4x7x1.6mm, smooth type, Myonic, spare part for dental handpieces, with steel cage and steel balls.

Handpiece compatibility

- Sirona T2 Line A40/L – T3 Line E40/L Contra-angle blue ring: push side
- Sirona T4 Line B40 Contra-angle blue ring: push side
- W&H WA-56 LT New Contra-angle blue ring: bur side and push side
- W&H Synea Fusion WG-56LT Old Model Contra-angle blue ring: push side
- W&H Synea Fusion WG-56LT New Model Contra-angle blue ring: bur side and push side
- W&H Synea Fusion WG-66 LT Old Model Contra-angle blue ring: push side

Technical features

- Dimensions: height 1.6mm, internal diameter 4mm, external diameter 7mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- Myonic brand
- 12-month warranty