

Bearing, low speed, 4x9x2,5 steel



Bearing for low speed handpiece, 4x9x2.5mm dimensions, smooth, steel cage, steel balls, Tecnomed Italia brand.

Shipping expected in approximately 15 days

SKU: CM414

Categories: Internal diameter 4-5

Description

Low speed handpiece bearing 4x9x2.5mm, smooth, by Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: 4x9x2.5mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Tecnomed Italia brand
- 12-month warranty