

## Bearing, low speed, 4x9x4 steel, ZZ,flanged



Bearing for low speed handpiece, 4x9x4mm dimensions, smooth flanged, steel cage, steel balls, closed on both sides, Tecnomed Italia brand.

**Available from 22/06/2026**

**SKU:** CM421F

**Categories:** Internal diameter 4-5

### Description

---

Low speed handpiece bearing 4x9x4mm, smooth flanged, by Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

### Technical features

- Dimensions: 4x9x4mm
- Smooth flanged bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty