


## Bearing, low speed, 5x11x3 steel

Bearing for low speed handpiece, 5x11x3mm dimensions, smooth, steel cage, steel balls, Tecnomed Italia brand.



**Available**

**SKU:** CM532

**Categories:** Internal diameter 4-5

### Description

---

Low speed handpiece bearing 5x11x3mm, smooth, by Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

### Technical features

- Dimensions: 5x11x3mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Tecnomed Italia brand
- 12-month warranty