


Bearing, low speed, 5x8x2 steel, Z,flanged



Bearing for low speed handpiece, 5x8x2mm dimensions, smooth flanged, steel cage, steel balls, closed on one side, Tecnomed Italia brand.

Shipping expected in approximately 15 days

SKU: CM523

Categories: Internal diameter 4-5

Description

Low speed handpiece bearing 5x8x2mm, smooth flanged, by Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: 5x8x2mm
- Smooth flanged bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Closed on one side
- Tecnomed Italia brand
- 12-month warranty