

Bearing, low speed, 5x8x2 steel



Bearing for low speed handpiece, dimensions 5x8x2mm, smooth, steel cage, steel balls, Myonic brand, compatible with KaVo, Ivory, Bien-Air and NSK.

Available

SKU: CM512

Categories: Internal diameter 4-5

Description

Low speed handpiece bearing, 5x8x2mm, smooth type, Myonic, spare part for dental handpieces, with steel cage and steel balls.

Handpiece compatibility

- KaVo SURGmatic S201XC/XL green ring: push side
- KaVo SURGmatic S201L green ring: push side
- Ivory SL4201C 20:1 green ring: push side
- Bien-Air CA 1:1 L / Diamline blue ring: push side
- Bien-Air CA 20:1 L / Diamline green ring: push side
- NSK Ti-Max X-SG20L/DSG20L green ring: push side
- NSK Ti-Max SG20 green ring: push side
- KaVo C3 – CL3 – Ivory SL631C Head Green Ring 3:1: push side

Technical features

- Dimensions: height 2mm, internal diameter 5mm, external diameter 8mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- Myonic brand
- 12-month warranty