

Bearing, low speed, 5x8x2,5 steel, ZZ, flanged



Bearing for low speed handpiece, 5x8x2.5mm dimensions, smooth flanged, steel cage, steel balls, closed on both sides, Tecnomed Italia brand.

Available

SKU: CM563M

Categories: Internal diameter 4-5

Description

Low speed handpiece bearing 5x8x2.5mm, smooth flanged, by Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: 5x8x2.5mm
- Smooth flanged bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty