

Bearing, low speed, 5x9,525x3,175 steel, ZZ



Bearing for low speed handpiece, 5x9.525x3.175mm dimensions, smooth, steel cage, steel balls, closed on both sides, Tecnomed Italia brand.

Shipping expected in approximately 15 days

SKU: CM556

Categories: Internal diameter 4-5

Description

Low speed handpiece bearing 5x9.525x3.175mm, smooth, Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: 5x9.525x3.175mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty