

## Low speed handpiece bearing 7.938×15.875×6mm, smooth, steel balls



Bearing for low speed handpiece micromotor, 7.938×15.875×6mm dimensions, smooth, steel cage, steel balls, closed on both sides, Tecnomed Italia brand, compatible with KaVo K9, K11 and K12.

**Available**

**SKU:** CU111K

**Categories:** Internal diameter 6-12

### Description

---

Low speed handpiece micromotor bearing 7.938×15.875×6mm, smooth by Tecnomed Italia, spare bearing for laboratory micromotors, with steel cage and steel balls.

### Laboratory micromotor compatibility

- KaVo K9 4930 Blue
- KaVo K11
- KaVo K12

### Technical features

- Dimensions: 7.938×15.875×6mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty