

Low speed handpiece bearing 7×14.1×7.5mm, groove, steel balls



Bearing for low speed handpiece micromotor, 7×14.1×7.5mm dimensions, groove, steel cage, steel balls, closed on both sides, Tecnomed Italia brand, compatible with KaVo K4-K5, K9 and K12.

Available

SKU: CM727

Categories: Internal diameter 6-12

Description

Low speed handpiece micromotor bearing 7×14.1×7.5mm, groove, by Tecnomed Italia, spare bearing for laboratory micromotors, with steel cage and steel balls.

Laboratory micromotor compatibility

- KaVo K4 – K5
- KaVo K9 4930 Blue
- KaVo K12

Technical features

- Dimensions: 7×14.1×7.5mm
- Groove bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty