

## Low speed handpiece bearing 8×14×4mm, smooth, steel balls



Bearing for low speed handpiece, 8×14×4mm dimensions, smooth, steel cage, steel balls, closed on both sides, Tecnomed Italia brand, Kavo K9 compatible.

**Available**

**SKU:** K220148

**Categories:** Internal diameter 6-12

### Description

---

Low speed handpiece bearing 8×14×4mm, smooth, Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

### Compatibility

- KaVo K9 (EWL 955 Grey-Brown)

### Technical features

- Dimensions: 8×14×4mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty