


Bearing, low speed, 8x16x4 steel

Bearing for low speed handpiece, 8x16x4mm dimensions, smooth, steel cage, steel balls, Tecnomed Italia brand.



Available

SKU: CM834

Categories: Internal diameter 6-12

Description

Low speed handpiece bearing 8x16x4mm, smooth, by Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: 8x16x4mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Tecnomed Italia brand
- 12-month warranty