



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/bowie-dick-test-pack-of-20/>

## Bowie Dick test, pack of 20

Bowie Dick test for detecting abnormalities during the cycle, 20-pack

**89,27 €**

**Available**

**SKU:** DEA4.01

**Categories:** Sterilization test



### Description

---

Test Bowie & Dick is a physical test for B –class autoclaves.

This test is used to assess the ability of steam penetration into porous bodies.

The presence of air inside the autoclave poses a threat to the safety of the sterilization process. This problem is solved through the disposable ready package Bowie & Dick as required by the current regulations.

Unlike traditional devices, the use of the ready package has significant advantages that ensure the test standards are maintained for each test carried out and allow you to immediately detect the removal of air during the sterilization process.

Bowie & Dick test is made in the empty chamber after the vacuum test and before the normal use of the autoclave.

- It allows you to immediately and accurately detect the presence of any anomalies during the Bowie Dick cycle
- The resulting response remains unchanged over time
- The external package is equipped with a toning indicator that signals, at a glance, whether the package was used or not.
- The materials which it is made with are fully recyclable and the toning chemical reagent is non-toxic
- Each package of Bowie & Dick is equipped with a complete card for the easy-to-fill data storage.
- Pack containing 20 pieces.