



Cabinet 500x480x900h with ceramic washbasin, top 1500mm with splash guard

Available on request

SKU: WASHT15SO

Categories: Sinks

Description

The **clinical 500x480x900h cabinet with ceramic sink** and 1500mm top with backsplash is engineered to support the workflow of modern practices, providing an organised, hygienic and highly functional workstation. Its modular structure enables a rational layout, improving efficiency and enhancing the patient's perception of the clinical environment.

STRUCTURE AND MATERIALS

Constructed from electro-galvanised steel in RAL 9010 with a 10/10 thickness, the cabinet ensures stability and long-term reliability. Fronts and handles feature an antibacterial epoxy coating, ideal for maintaining hygiene in clinical settings. The flush-mounted ceramic sink integrates seamlessly with the worktop, ensuring practicality and easy sanitisation.

CLINICAL WORK SURFACES

The scratch-resistant and water-repellent laminated top provides a continuous and durable work surface developed for intensive clinical routines. Its generous width supports efficient instrument organisation and smooth workflow during procedures.

INTEGRATED CLINICAL EQUIPMENT

The sink module includes a **clinical lever mixer** for precise and safe temperature control, while the **electric foot control** allows hands-free water delivery, ensuring optimal hygiene. Anodised aluminium handles — available in eight colours — allow the cabinet to match the aesthetic identity of the practice.

DESIGN, FUNCTIONALITY AND ORGANISATION



Clean lines, compact geometry and uniform surfaces contribute to a modern and professional look. The modular layout ensures that all instruments remain well arranged and immediately accessible, supporting a smoother clinical workflow and improving operator comfort.