



## Cabinet 500x480x900h with ceramic washbasin, top 2000mm with splash guard

Available on request

SKU: WASHT20SO

Categories: Sinks

### Description

---

The **clinical cabinet with ceramic sink and 2000mm top** is designed to provide a highly organised, hygienic and efficient clinical workstation. Its modular configuration supports a rational layout, optimising workflow and enhancing the patient experience.

#### STRUCTURE AND MATERIALS

Manufactured from electro-galvanised steel in RAL 9010 with 10/10 thickness, the cabinet ensures solidity, stability and long-term durability. Fronts and handles coated with antibacterial epoxy powder help maintain hygiene in demanding clinical environments.

The flush-mounted ceramic sink offers seamless integration, allowing for easier cleaning and a continuous, practical working surface.

#### CLINICAL WORK SURFACE

The **2000mm** scratch-resistant and water-repellent laminated top provides an extended, durable and hygienic work area. Its generous length enhances instrument organisation, reduces operator movements and supports ergonomic clinical routines.

#### INTEGRATED CLINICAL EQUIPMENT

The sink module includes a **clinical lever mixer** for precise and safe control of water temperature, while the **electric foot control** enables hands-free operation, ensuring optimal hygiene. Anodised aluminium handles, available in eight colours, allow aesthetic customisation suited to the identity of the practice.

#### PROFESSIONAL DESIGN AND ORGANISATION



Clean lines, compact geometry and uniform surfaces create a modern and professional environment. Internal module organisation ensures that all instruments and materials remain accessible, supporting smooth and efficient clinical procedures.