




Cabinet 500x480x900h with ceramic washbasin, top 2000mm with splash guard



Available on request

SKU: WASHT20SO

Categories: Sinks

Description

The **clinical cabinet with ceramic sink and 2000mm top** is designed to provide a highly organised, hygienic and efficient clinical workstation. Its modular configuration supports a rational layout, optimising workflow and enhancing the patient experience.

STRUCTURE AND MATERIALS

Manufactured from electro-galvanised steel in RAL 9010 with 10/10 thickness, the cabinet ensures solidity, stability and long-term durability. Fronts and handles coated with antibacterial epoxy powder help maintain hygiene in demanding clinical environments.

The flush-mounted ceramic sink offers seamless integration, allowing for easier cleaning and a continuous, practical working surface.

CLINICAL WORK SURFACE

The **2000mm** scratch-resistant and water-repellent laminated top provides an extended, durable and hygienic work area. Its generous length enhances instrument organisation, reduces operator movements and supports ergonomic clinical routines.

INTEGRATED CLINICAL EQUIPMENT

The sink module includes a **clinical lever mixer** for precise and safe control of water temperature, while the **electric foot control** enables hands-free operation, ensuring optimal hygiene. Anodised aluminium handles, available in eight colours, allow aesthetic customisation suited to the identity of the practice.

PROFESSIONAL DESIGN AND ORGANISATION



Clean lines, compact geometry and uniform surfaces create a modern and professional environment. Internal module organisation ensures that all instruments and materials remain accessible, supporting smooth and efficient clinical procedures.