




Printed on 05/06/2025

<https://shop.tecnomeditalia.com/en/product/cabinet-for-practice-glass-top-50cm-standard-with-left-side-sink/>

Cabinet, for practice, glass top 50cm, standard, with left side sink



Cabinet for medical practice, glass top 50cm, made of robust electro-galvanised sheet steel, complete with left side sink.

5.701,70 €

Available on request

SKU: DE9.21SXV

Categories: Studio

Description

Dental cabinetry left side sink glass.

Our range features top quality materials and simple designs for maximum versatility.

Tecnomed cabinets are superbly ergonomic and designed to make the lives of working dentists easier.

Made from robust, electro-galvanised sheet steel in a uniform, 10/10 thickness, these cabinets are quick and easy to clean and designed to stop dust and debris building up in hollows and cracks.

Technical features Dental cabinetry leftside sink glass DE9.21SXV.

Top	Glass
Washbasin Cabinet	•
4 compartment drawer unit	•
4 compartment service unit	•
Full door unit	• 1 glass shelf
Washbasin – material	Left – Glass
Lever type medical tap	•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	No
Warranty	24 mesi
CE	•