



Cabinet, for practice, glass top 50cm, standard, with left side sink

Cabinet for medical practice, glass top 50cm, made of robust electro-galvanised sheet steel, complete with left side sink.

Available on request

SKU: DE9.21SXV

Categories: Cabinet Compositions

Description

The **clinic cabinet with glass top and left-side sink** is designed for dental practices, medical clinics, veterinary rooms and beauty studios.

Made of homogeneous electro-galvanised steel (10/10 thickness), it ensures durability, structural stability and easy cleaning in demanding clinical environments.

Ergonomic and functional design

Its compact and linear layout optimises daily workflow, supporting smooth and comfortable movements for operators.

The left-side sink ensures natural and immediate access during routine procedures.

High-quality hygienic materials

The cabinet features a RAL 9010 finish, while the fronts and handles are coated with antibacterial epoxy powders to minimise dust accumulation and simplify sanitisation.

The **scratch-resistant glass top** provides a durable, hygienic and easy-to-clean surface, perfectly suited to intensive clinical use.

Custom configuration – Made in Italy

Each cabinet can be customised to match room layout and operational needs, ensuring versatility and optimal space management.

Entirely Made in Italy, it reflects meticulous craftsmanship, reliable performance and long-term durability.

This clinic cabinet with glass top and left-side sink is the ideal furnishing solution for modern healthcare environments requiring hygiene, organisation and ergonomic design.



Technical features.

Top	Glass
Washbasin Cabinet	<ul style="list-style-type: none">•
4 compartment drawer unit	<ul style="list-style-type: none">•
4 compartment service unit	<ul style="list-style-type: none">•
Full door unit	<ul style="list-style-type: none">• 1 glass shelf
Washbasin – material	Left – Glass
Lever type medical tap	<ul style="list-style-type: none">•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	No
Warranty	24 mesi
CE	<ul style="list-style-type: none">•