


Cabinet, for practice, glass top 50cm, standard, with right side sink



Cabinet for medical practice, glass top 50cm, made of robust electro-galvanised sheet steel, complete with right side sink.

Available on request

SKU: DE9.21DXV

Categories: Cabinet Compositions

Description

The **clinic cabinet with glass top and right-side sink** is designed for dental practices, medical clinics, veterinary rooms and beauty studios. Made of homogeneous electro-galvanised steel (10/10 thickness), it ensures strength, long-term durability and easy cleaning in demanding clinical environments.

Ergonomic and functional design

Its compact and linear structure optimises daily workflow and supports efficient, comfortable movements for operators. The right-side sink provides natural and immediate access during routine procedures.

High-quality hygienic materials

The cabinet features a RAL 9010 finish, with fronts and handles coated in antibacterial epoxy powder to reduce dust build-up and simplify sanitisation. The scratch-resistant glass top enhances hygiene, offering a resistant, easy-to-clean surface suitable for clinical settings.

Custom configuration – Made in Italy

Each unit can be customised to match room layout and operational requirements, ensuring flexibility and optimal use of space. Completely Made in Italy, it reflects high-quality craftsmanship, attention to detail and long-lasting reliability.

This clinic cabinet with glass top and right-side sink is the perfect furnishing solution for modern healthcare environments requiring hygiene, organisation and ergonomic design.

Technical features Dental cabinetry rightside sink.

Top	Glass
-----	-------



Washbasin Cabinet	<ul style="list-style-type: none">•
4 compartment drawer unit	<ul style="list-style-type: none">•
4 compartment service unit	<ul style="list-style-type: none">•
Full door unit	<ul style="list-style-type: none">• 1 glass shelf
Washbasin – material	Right – Glass
Lever type medical tap	<ul style="list-style-type: none">•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	No
Warranty	24 mesi
CE	<ul style="list-style-type: none">•