




Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/cabinet-for-practice-glass-top-50cm-standard-with-right-side-sink/>

Cabinet, for practice, glass top 50cm, standard, with right side sink



Cabinet for medical practice, glass top 50cm, made of robust electro-galvanised sheet steel, complete with right side sink.

5.701,70 €

Available on request

SKU: DE9.21DXV

Categories: Studio

Description

Dental cabinetry rightside sink glass.

Our range features top quality materials and simple designs for maximum versatility.

Tecnomed cabinets are superbly ergonomic and designed to make the lives of working dentists easier.

Made from robust, electro-galvanised sheet steel in a uniform, 10/10 thickness, these cabinets are quick and easy to clean and designed to stop dust and debris building up in hollows and cracks.

Technical features Dental cabinetry rightside sink + DE9.21DXV.

Top	Glass
Washbasin Cabinet	•
4 compartment drawer unit	•
4 compartment service unit	•
Full door unit	• 1 glass shelf
Washbasin – material	Right – Glass
Lever type medical tap	•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	No
Warranty	24 mesi
CE	•