



Cabinet for practice, top 50 cm, standard, complete with left side sink

Cabinet for medical practice, top 50cm, made of robust electrogalvanised sheet steel, complete with left side sink.

Available on request

SKU: DE9.21SX

Categories: Cabinet Compositions

Description

The **Dental Cabinetry with Left-Side Sink DE9.21SX** combines elegant design, premium materials, and ergonomic functionality to support the daily needs of dental professionals.

Manufactured from robust electro-galvanised sheet steel with a uniform 10/10 thickness, this cabinetry ensures durability, hygiene, and long-lasting performance. Its smooth surfaces are quick and easy to clean, while the design eliminates hollows and cracks that could accumulate dust or debris.

Tecnomed cabinetry is engineered for maximum versatility and comfort. With a focus on ergonomics, these cabinets make workflows more efficient, creating a safe, practical, and organized environment in the dental clinic.

Technical features Dental cabinetry left-side sink DE9.21SX.

Тор	Laminate
Washbasin Cabinet	•
4 compartment drawer unit	•
4 compartment service unit	•
Full door unit	• 1 glass shelf
Washbasin – material	Left – ABS



Lever type medical tap	•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	No
Warranty	24 mesi
CE	•