



Cabinet for practice, top 50 cm, standard, complete with left side sink

Cabinet for medical practice, top 50cm, made of robust electro-galvanised sheet steel, complete with left side sink.

Available on request

SKU: DE9.21SX

Categories: Cabinet Compositions

Description

The **clinic cabinet with left-side sink**, equipped with a 50 cm top, is designed for **dental practices, medical clinics, veterinary rooms and beauty studios**.

Made of **homogeneous electro-galvanised steel** (10/10 thickness), it ensures strength, durability and easy cleaning in high-demand clinical environments.

Ergonomic and functional design

Its **compact and linear design** enhances workflow, allowing operators to work comfortably and efficiently. The right-side sink provides immediate, natural access during everyday procedures.

High-quality hygienic materials

The structure features a **RAL 9010 finish**, while the fronts and handles are coated with **antibacterial epoxy powders**, reducing dust accumulation and simplifying daily sanitisation.

Custom configuration – Made in Italy

The cabinet can be **customised** to suit any room layout, offering versatility and adaptability. Completely **Made in Italy**, it embodies high-quality manufacturing, attention to detail and long-term reliability.

This **clinic cabinet with left-side sink** is the perfect furnishing solution for modern healthcare environments requiring hygiene, organisation and ergonomic design.

Technical features

Top	Laminate
-----	----------



Washbasin Cabinet	<ul style="list-style-type: none">•
4 compartment drawer unit	<ul style="list-style-type: none">•
4 compartment service unit	<ul style="list-style-type: none">•
Full door unit	<ul style="list-style-type: none">• 1 glass shelf
Washbasin – material	Left – ABS
Lever type medical tap	<ul style="list-style-type: none">•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	No
Warranty	24 mesi
CE	<ul style="list-style-type: none">•