



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/cabinet-for-practice-top-50-cmleftside-sink-solenoid-valve/>



Cabinet for practice, top 50 cm, leftside sink + solenoid valve

Cabinet for medical practice with solenoid valve, top 50cm, made of robust electro-galvanised sheet steel, complete with left side sink.

3.043,04 €

Available on request

SKU: DE9.24SX

Categories: Studio

Description

Dental cabinetry leftside sink + solenoid valve.

Our range features top quality materials and simple designs for maximum versatility.

Tecnomed cabinets are superbly ergonomic and designed to make the lives of working dentists easier.

Made from robust, electro-galvanised sheet steel in a uniform, 10/10 thickness, these cabinets are quick and easy to clean and designed to stop dust and debris building up in hollows and cracks.

Technical features Dental cabinetry leftside sink + solenoid valve DE9.24SX.

Top	Laminate
Washbasin Cabinet	•
4 compartment drawer unit	•
4 compartment service unit	•
Full door unit	• 1 glass shelf
Washbasin – material	Left – ABS
Lever type medical tap	•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	•
Warranty	24 mesi
CE	•