



## Cabinet for practice, top 50 cm, left side sink + solenoid valve

Cabinet for medical practice with solenoid valve, top 50cm, made of robust electro-galvanised sheet steel, complete with left side sink.

**Available on request**

**SKU:** DE9.24SX

**Categories:** Cabinet Compositions

### Description

---

The **clinic cabinet with solenoid valve and left-side sink**, equipped with a 50 cm top, is designed for **dental practices, medical clinics, veterinary rooms and beauty studios**.

Made of **homogeneous electro-galvanised steel** (10/10 thickness), it ensures strength, durability and easy cleaning in high-demand clinical environments.

### Ergonomic and functional design

Its **compact and linear design** enhances workflow, allowing operators to work comfortably and efficiently. The right-side sink provides immediate, natural access during everyday procedures.

### High-quality hygienic materials

The structure features a **RAL 9010 finish**, while the fronts and handles are coated with **antibacterial epoxy powders**, reducing dust accumulation and simplifying daily sanitisation.

The integrated **solenoid valve** guarantees precise water control, improving hygiene and ensuring safe operation in clinical settings.

### Custom configuration – Made in Italy

The cabinet can be **customised** to suit any room layout, offering versatility and adaptability. Completely **Made in Italy**, it embodies high-quality manufacturing, attention to detail and long-term reliability.

This **clinic cabinet with solenoid valve and left-side sink** is the perfect furnishing solution for modern healthcare environments requiring hygiene, organisation and ergonomic design.

### Technical features



Top	Laminate
Washbasin Cabinet	<ul style="list-style-type: none"><li>•</li></ul>
4 compartment drawer unit	<ul style="list-style-type: none"><li>•</li></ul>
4 compartment service unit	<ul style="list-style-type: none"><li>•</li></ul>
Full door unit	<ul style="list-style-type: none"><li>• 1 glass shelf</li></ul>
Washbasin – material	Left – ABS
Lever type medical tap	<ul style="list-style-type: none"><li>•</li></ul>
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	<ul style="list-style-type: none"><li>•</li></ul>
Warranty	24 mesi
CE	<ul style="list-style-type: none"><li>•</li></ul>