



Cabinet for practice, top 50 cm, right side sink + solenoid valve

Cabinet for medical practice with solenoid valve, top 50cm, made of robust electro-galvanised sheet steel, complete with right side sink.

Available on request

SKU: DE9.24DX

Categories: Cabinet Compositions

Description

The **clinic cabinet with solenoid valve and right-side sink**, equipped with a 50 cm top, is designed for **dental practices, medical clinics, veterinary rooms and beauty studios**.

Made of **homogeneous electro-galvanised steel** (10/10 thickness), it ensures strength, durability and easy cleaning in high-demand clinical environments.

Ergonomic and functional design

Its **compact and linear design** enhances workflow, allowing operators to work comfortably and efficiently. The right-side sink provides immediate, natural access during everyday procedures.

High-quality hygienic materials

The structure features a **RAL 9010 finish**, while the fronts and handles are coated with **antibacterial epoxy powders**, reducing dust accumulation and simplifying daily sanitisation.

The integrated **solenoid valve** guarantees precise water control, improving hygiene and ensuring safe operation in clinical settings.

Custom configuration – Made in Italy

The cabinet can be **customised** to suit any room layout, offering versatility and adaptability.

Completely **Made in Italy**, it embodies high-quality manufacturing, attention to detail and long-term reliability.

This **clinic cabinet with solenoid valve and right-side sink** is the perfect furnishing solution for modern healthcare environments requiring hygiene, organisation and ergonomic design.

Technical features.



Top	Laminate
Washbasin Cabinet	•
4 compartment drawer unit	•
4 compartment service unit	•
Full door unit	• 1 glass shelf
Washbasin – material	Right – ABS
Lever type medical tap	•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	•
Warranty	24 mesi
CE	•