



## Cabinet, for practice, top 50cm, standard, complete with right side sink

Cabinet for medical practice, top 50cm, made of robust electro-galvanised sheet steel, standard, complete with right side sink.

**Estimated shipping in 30 days**

**SKU:** DE9.21DX

**Categories:** Cabinet Compositions

### Description

The **clinic cabinet with right-side sink**, equipped with a 50 cm top, is designed for **dental practices, medical clinics, veterinary rooms and beauty studios**.

Made of **homogeneous electro-galvanised steel** (10/10 thickness), it ensures strength, durability and easy cleaning in high-demand clinical environments.

### Ergonomic and functional design

Its **compact and linear design** enhances workflow, allowing operators to work comfortably and efficiently. The right-side sink provides immediate, natural access during everyday procedures.

### High-quality hygienic materials

The structure features a **RAL 9010 finish**, while the fronts and handles are coated with **antibacterial epoxy powders**, reducing dust accumulation and simplifying daily sanitisation.

### Custom configuration – Made in Italy

The cabinet can be **customised** to suit any room layout, offering versatility and adaptability. Completely **Made in Italy**, it embodies high-quality manufacturing, attention to detail and long-term reliability.

This **clinic cabinet with right-side sink** is the perfect furnishing solution for modern healthcare environments requiring hygiene, organisation and ergonomic design.

Technical features.

Top	Laminate
-----	----------



Washbasin Cabinet	<ul style="list-style-type: none"><li>•</li></ul>
4 compartment drawer unit	<ul style="list-style-type: none"><li>•</li></ul>
4 compartment service unit	<ul style="list-style-type: none"><li>•</li></ul>
Full door unit	<ul style="list-style-type: none"><li>• 1 glass shelf</li></ul>
Washbasin – material	Right – ABS
Lever type medical tap	<ul style="list-style-type: none"><li>•</li></ul>
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Electrovalve	No
Warranty	24 mesi
CE	<ul style="list-style-type: none"><li>•</li></ul>