



## Cabinet, for practice, top 50cm, standard, complete with right side sink

Cabinet for medical practice, top 50cm, made of robust electro-galvanised sheet steel, standard, complete with right side sink.

**Available on request**

**SKU:** DE9.21DX

**Categories:** Studio

### Description

Dental cabinetry rightside sink.

Our range features top quality materials and simple designs for maximum versatility.

Tecnomed cabinets are superbly ergonomic and designed to make the lives of working dentists easier.

Made from robust, electro-galvanised sheet steel in a uniform, 10/10 thickness, these cabinets are quick and easy to clean and designed to stop dust and debris building up in hollows and cracks.

Technical features Dental cabinetry rightside sink + solenoid valve DE9.21DX.

Top	Laminate
Washbasin Cabinet	•
4 compartment drawer unit	•
4 compartment service unit	•
Full door unit	• 1 glass shelf
Washbasin – material	Right – ABS
Lever type medical tap	•
Worktop dimensions	2080 x 600 mm



Height from floor	900 mm
Electrovalve	No
Warranty	24 mesi
CE	•