



Cabinet sterilization top 60 cm, complete with leftside sink

Sterilization cabinet, top 60cm, made of robust electro-galvanised sheet steel, complete with left side sink.

Available on request

SKU: DE9.21SX-60

Categories: Sterilization Cabinets

Description

The **sterilisation cabinet with left-hand sink** is designed to deliver **versatility, hygiene and functionality** in clinical and dental settings. Every feature reflects an ergonomic concept that enhances workflow efficiency and supports comfortable, precise movements throughout daily operations.

Built from **electro-galvanised steel 10/10 gauge**, the cabinet provides uniform strength and smooth, seamless surfaces that minimise dust and residue accumulation. This ensures **maximum cleanliness** and simplifies maintenance, helping keep clinical environments hygienic and professional.

The **laminate worktop** includes a **left-hand ABS sink**, offering a spacious, durable surface for washing and disinfecting instruments. The layout features an **under-sink unit**, **two four-compartment drawer modules**, and a **cabinet with a hinged door and glass shelf**, ensuring practical organisation and quick access to instruments and materials during sterilisation cycles.

At a height of **900 mm**, the cabinet promotes ergonomic posture and effortless handling. The product is **CE-certified** and covered by a **24-month warranty**, representing a balanced combination of **durability, hygiene and ergonomic performance**—ideal for sterilisation areas, dental practices, clinics and laboratories requiring professional-grade, easy-to-clean furniture.

Technical features.

Top	Laminate
Washbasin Cabinet	•
4 compartment drawer unit	• x2



Full door unit	<ul style="list-style-type: none">• 1 glass shelf
Washbasin – material	left – ABS
Lever type medical tap	<ul style="list-style-type: none">•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Warranty	24 mesi
CE	<ul style="list-style-type: none">•