



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/cabinet-sterilization-top-60-cm-complete-with-leftside-sink/>



Cabinet sterilization top 60 cm, complete with leftside sink

Sterilization cabinet, top 60cm, made of robust electro-galvanised sheet steel, complete with left side sink.

2.818,82 €

Available on request

SKU: DE9.21SX-60

Categories: Sterilization

Description

Sterilization cabinet leftside sink.

Our range features top quality materials and simple designs for maximum versatility. Tecnomed cabinets are superbly ergonomic and designed to make the lives of working dentists easier. Made from robust, electro-galvanised sheet steel in a uniform, 10/10 thickness, these cabinets are quick and easy to clean and designed to stop dust and debris building up in hollows and cracks.

Technical features Sterilization cabinet leftside sink.

Top	Laminate
Washbasin Cabinet	•
4 compartment drawer unit	• x2
Full door unit	• 1 glass shelf
Washbasin – material	left – ABS
Lever type medical tap	•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Warranty	24 mesi
CE	•