


Cabinet sterilization top 60 cm, complete with rightside sink



Sterilization cabinet, top 60cm, made of robust electro-galvanised sheet steel, complete with right side sink.

Available on request

SKU: DE9.21DX-60

Categories: Sterilization Cabinets

Description

The **sterilisation cabinet with right-hand sink** emerges from a design line developed for **versatility, hygiene and ergonomics** in clinical and dental environments. The structure is inspired by a concept of high quality, engineered to enhance operator movement and promote the true essence of ergonomic clinical furnishings.

Constructed in **electro-galvanised steel 10/10 gauge**, the cabinet ensures uniform strength and durability while providing smooth surfaces that inhibit dust and residue accumulation in micro-indentations. This design fosters **maximum hygiene** and facilitates quicker, more efficient cleaning and maintenance.

The **laminate worktop** integrates a **right-hand sink in ABS**, providing a spacious and durable surface ideal for instrument washing and disinfection. The layout includes an **under-sink cupboard**, **two four-compartment drawer units**, and a **cabinet with hinged door and glass shelf**, offering well-organised storage for tools and materials and streamlining workflow during sterilisation tasks.

With a height of **900 mm**, the unit supports comfortable posture and easy access to storage compartments. The product is **CE-certified** and covered by a **24-month warranty**. It represents an ideal balance of **durability, ergonomics and cleanliness**, suited for sterilisation rooms, dental practices, clinics and laboratories requiring professional-grade furniture.

Technical features Sterilization cabinet rightside sink.

Top	Laminate
Washbasin Cabinet	•



4 compartment drawer unit	<ul style="list-style-type: none">• x2
Full door unit	<ul style="list-style-type: none">• 1 glass shelf
Washbasin – material	Right – ABS
Lever type medical tap	<ul style="list-style-type: none">•
Worktop dimensions	2080 x 600 mm
Height from floor	900 mm
Warranty	24 mesi
CE	<ul style="list-style-type: none">•