



Easy Tray inox surgical equipment trolley, 5 removable trays

Easy Tray white surgical equipment trolley, with 5 stainless steel removable and sterilizable trays, dimensions 46,5x31,5x2 cm.

Available

SKU: OF1003BI.00

Categories: Equipment Carts

Description

Easy Tray White is a powered medical trolley manufactured in white powder-coated steel, designed for the organised placement and power supply of electromedical, dental and surgical equipment in high-intensity clinical and operating environments. Finished in white powder coating, Easy Tray White blends naturally into modern healthcare environments while delivering the same functional performance as the stainless steel version.

Designed for dental practices, outpatient clinics, operating theatres and hospital environments, the 5-tray configuration with integrated height-adjustable IV pole provides a fully equipped mobile workstation capable of supporting multiple devices and instruments during clinical procedures.

With removable, sterilizable trays, instrument management becomes clear and intuitive, improving ergonomics and workflow efficiency. Thanks to the integrated power system, equipment can be powered directly on the trolley, reducing cable clutter and increasing operational safety.

Concealed electrical cable management keeps the work area clean, orderly and easy to disinfect. The **solid and stable structure** makes Easy Tray White suitable for intensive clinical use, ensuring long-term reliability without unnecessary complexity.

Easy Tray White cart is also a **modular platform**, compatible with a wide range of optional accessories, allowing configuration according to the specific operational needs of the professional user.

Key features

- Powered medical trolley for clinical and electromedical equipment
- White powder-coated steel structure, robust and easy to clean and disinfect
- **5 removable and sterilizable trays** for instruments and devices (465 x 315 x 20 mm)



- **Height-adjustable IV pole included**
- Concealed electrical cable management for improved safety and organisation
- Strong, stable design suitable for intensive clinical use
- Medical-grade castors with brakes for stability during procedures
- Suitable as a medical trolley for clinics, dental practices and surgical environments

Overall dimensions:

Height 145cm (215cm with IV pole fully extended) × Width 60cm × Depth 50cm

Optional accessories (available separately)

- Monitor support with VESA 75/100 mount for medical monitors (ref. SVET9010L)
- Handpiece support arm (ref. DE1003X.5M)
- Implantology handpiece holder (ref. DE1003X.1)
- Piezosurgery handpiece holder (ref. DE1003X.2)
- Electrosurgical handpiece holder (ref. DE1003X.3)
- Medical camera arm for OKKIO camera system (ref. DE1003X.6K)
- Additional tray 415 × 300 × 20 mm (ref. DE1003X.02)