



Easy Tray inox surgical equipment trolley, 5 removable trays

Easy Tray Inox surgical equipment trolley, made of stainless steel with 5 removable and sterilizable trays, dimensions 46,5x31,5x2 cm.

Available

SKU: OF1003X.00

Categories: Equipment Carts

Description

Easy Tray Inox surgical equipment trolley **is the best solution for host you medical equipment in a orderly manner.**

Versatile and modular: Made entirely of electro galvanized white enamelled sheet, this cart is modular and configurable with different accessories according to your needs.

No cable in your way: you'll get a cleaner, more sterile environment thanks to hidden electrical connections and a concealed power supply.

One cart for all your equipment: 5 sterilizable and removable trays for host complementary surgical equipment such as scalers, electrosurgery units, piezoelectric surgical devices, implant micromotors, endodontic micromotors, lasers, surgical Stryker equipment, whitening devices, multi-parameter monitors. A single, integrated, height-adjustable infusion set stand can hold up to three bags of saline solution to ensure easy operation.

Fields of use: ideal for medical, surgical and dental settings, providing a safe and well-organised way to manage clinical equipment.

Technical feature Easy Tray Inox surgical equipment trolley:

- 5 stainless steel tray (46,5x31,5x2 cm), removable and sterilizable;
- height-adjustable infusion set stand
- Centralized power system for equipment.