



Cart Easy Tray inox with 5 tray

Available on request

SKU: OF1003X.00

Categories: Equipment Carts

Description

Easy Tray Inox 5 trays cart for equipment is the best solution for host you medical equipment in a orderly manner.

Versatile and modular: Made entirely of electro galvanized white enamelled sheet, this cart is modular and configurable with different accessories according to your needs.

No cable in your way: you'll get a cleaner, more sterile environment thanks to hidden electrical connections and a concealed power supply.

One cart for all your equipment: 5 sterilizable and removable trays for host complementary surgical equipment such as scalers, electrosurgery units, piezoelectric surgical devices, implant micromotors, endodontic micromotors, lasers, surgical Stryker equipment, whitening devices, multi-parameter monitors. A single, integrated, height-adjustable infusion set stand can hold up to three bags of saline solution to ensure easy operation.

Technical feature Easy Tray Inox equipment cart:

- 5 stainless steel tray (46,5×31,5×2 cm), removable and sterilizable;
- height-adjustable infusion set stand
- Centralized power system for equipment.

