



## Easy Tray inox surgical equipment trolley, 5 removable trays

Easy Tray Inox surgical equipment trolley, made of stainless steel with 5 removable and sterilizable trays, dimensions 46,5x31,5x2 cm.

**Available**

**SKU:** OF1003X.00

**Categories:** Equipment Carts

### Description

---

**Easy Tray Inox** is a powered medical trolley manufactured in stainless steel, designed for the organised placement and power supply of electromedical, dental and surgical equipment in high-intensity clinical and operating environments.

The **5-tray configuration with integrated adjustable IV pole** is intended for dental practices, outpatient clinics, operating theatres and hospital settings where a fully equipped mobile workstation is required, capable of supporting multiple devices and instruments during clinical procedures.

The **removable, sterilizable trays** allow clear and intuitive instrument management, improving ergonomics and workflow efficiency. Thanks to the **integrated power system**, equipment can be powered directly on the trolley, reducing cable clutter and increasing operational safety.

Hidden electrical cable management keeps the work area clean, orderly and easy to disinfect. The **solid and stable stainless steel structure** makes Easy Tray Inox suitable for intensive clinical use, ensuring long-term reliability without unnecessary complexity.

Easy Tray Inox cart is also a **modular platform**, compatible with a wide range of optional accessories, allowing configuration according to the specific operational needs of the professional user.

---

### Key features

- Powered medical trolley for clinical and electromedical equipment
- Stainless steel structure, robust and easy to clean and disinfect



- **5 removable and sterilizable trays** for instruments and devices (465 x 315 x 20 mm)
- **Height-adjustable IV pole included**
- Concealed electrical cable management for improved safety and organisation
- Strong, stable design suitable for intensive clinical use
- Medical-grade castors with brakes for stability during procedures
- Suitable as a medical trolley for clinics, dental practices and surgical environments

**Overall dimensions:**

Height 145cm (215cm with IV pole fully extended) x Width 60cm x Depth 50cm

---

**Optional accessories (available separately)**

- Monitor support with VESA 75/100 mount for medical monitors (ref. SVET9010L)
- Handpiece support arm (ref. DE1003X.5M)
- Implantology handpiece holder (ref. DE1003X.1)
- Piezosurgery handpiece holder (ref. DE1003X.2)
- Electrosurgical handpiece holder (ref. DE1003X.3)
- Additional tray 415 x 300 x 20 mm (ref. DE1003X.02)