



Cart extension for Easy Tray Nano with 2 stainless trays 465x315mm

Extension for Easy Tray Nano trolley with two 465x315mm stainless steel trays to convert it into a 5-tray Easy Tray configuration.

Available

SKU: OF1003BIN.02

Categories: Equipment Carts

Description

The **Easy Tray Nano cart extension** with 2 stainless steel trays 465x315mm is designed to transform the Easy Tray Nano trolley (code OF1003BIN.01) into a 5-tray Easy Tray configuration.

The extension increases the operational capacity of the workstation while maintaining order, continuity of use and structural consistency with the original trolley.

Once configured as a 5-tray Easy Tray, the trolley can be further customised with optional accessories to suit specific clinical requirements.

Key features Cart extension Easy Tray Nano

- Extension compatible with Easy Tray Nano trolley (code OF1003BIN.01)
- Converts the trolley into a 5-tray Easy Tray configuration
- Includes 2 stainless steel instrument trays
- Designed for optional accessory integration
- Easy-to-clean surfaces
- Tray dimensions: 465x315mm