



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/conservative-scaler-tip-ems-compatible-rt1-woodpecker-p4d/>



## Conservative scaler tip, EMS compatible RT1, Woodpecker P4D

Diamond disc conservative drill for canal localisation and coronoradicular preparation. Compatib. EMS RT1, Woodpecker P4D

**44,00 €**

**Available**

**SKU:** PG053

**Categories:** Conservative and prosthetic

### Description

---

Conservative scaler tip PG053.

The range Conservative & Prosthetics of Tecnomed Italia inserts is ideal for the microinvasive preparation of cavities and for decementation.

Our diamond conservative tips let you target lesions with great accuracy without damaging healthy tissue.

They are also excellent for removing carious dentine without damaging healthy adjacent teeth.

Our inserts for the decementation of implants and bridges are specially designed to help you work more rapidly and effectively.

Diamond disc conservative drill for canal localisation and coronoradicular preparation.

Compatibility Conservative scaler tip PG053:

- EMS RT1;
- Woodpecker P4D.