



Low speed handpiece bearing 2x5x1.5mm, smooth, steel balls, 200000 rpm

Bearing for low speed handpiece, dimensions 2x5x1.5mm, smooth, steel cage, steel balls, closed on both sides, 200.000 rpm, Tecnomed Italia brand.

Available

SKU: CM2834

Categories: Internal diameter 2-3

Description

Low speed handpiece bearing, 2x5x1.5mm, smooth type, Tecnomed Italia, spare part for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: height 1.5mm, internal diameter 2mm, external diameter 5mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- 200.000 rpm
- Tecnomed Italia brand
- 12-month warranty