

## Cuscinetto bassa velocità 4x11x4 sszz li tmi gabbia sintetica rossa



Bearing for low speed handpiece, 4x11x4mm dimensions, smooth, steel cage, steel balls, closed on both sides, Tecnomed Italia brand.

**Shipping expected in approximately 15 days**

**SKU:** CM0174

**Categories:** Internal diameter 4-5

### Description

---

Low speed handpiece bearing 4x11x4mm, smooth, Tecnomed Italia, spare bearing for dental handpieces, with steel cage and steel balls.

### Technical features

- Dimensions: 4x11x4mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty