



Equipment trolley, white, 3 shelves, 600x410x890 mm (lpxh)

The ES1003BS Dental Cart is a fundamental element for any medical practice or laboratory thanks to its versatility and user-friendliness.

600X410X890 MM, white, 3 layers

Available

SKU: ES1003BS

Categories: Carts

Description

White Soft Trolley for Equipment 3 Shelves – Durable and Versatile Design for Clinics and Studios

The **white soft trolley for equipment 3 shelves** is an essential accessory for medical, dental and aesthetic clinics. Designed for strength and practicality, it provides a stable and reliable surface for instruments and devices in daily professional use.

Strong and durable structure.

Made of high-resistance metal, it ensures long-lasting stability even under frequent use. Its smooth surfaces are easy to clean and disinfect, maintaining strict hygiene standards in clinical environments.

Efficient storage and organisation.

Featuring **3 wide shelves**, this trolley allows orderly arrangement of equipment and tools. The **middle shelf** is height-adjustable with a secure locking system, ensuring maximum flexibility for different working setups.

Smooth mobility and clean design.

The castors provide effortless, quiet movement, allowing easy repositioning within tight spaces. The soft white finish adds a clean, professional look that complements any clinical interior.

Dimensions:

Trolley 60x41x89cm

Top shelf 53x29.5cm – middle 27.5x47.5cm – bottom 54x34cm