



## Equipment trolley, white, 3 shelves, 600x410x890 mm (lpxh)

The ES1003BS Dental Cart is a fundamental element for any medical practice or laboratory thanks to its versatility and user-friendliness.

600X410X890 MM, white, 3 layers

**Available**

**SKU:** ES1003BS

**Categories:** Carts

### Description

---

#### **White Soft Trolley for Equipment 3 Shelves – Durable and Versatile Design for Clinics and Studios**

The **white soft trolley for equipment 3 shelves** is an essential accessory for medical, dental and aesthetic clinics. Designed for strength and practicality, it provides a stable and reliable surface for instruments and devices in daily professional use.

##### **Strong and durable structure.**

Made of high-resistance metal, it ensures long-lasting stability even under frequent use. Its smooth surfaces are easy to clean and disinfect, maintaining strict hygiene standards in clinical environments.

##### **Efficient storage and organisation.**

Featuring **3 wide shelves**, this trolley allows orderly arrangement of equipment and tools. The **middle shelf** is height-adjustable with a secure locking system, ensuring maximum flexibility for different working setups.

##### **Smooth mobility and clean design.**

The castors provide effortless, quiet movement, allowing easy repositioning within tight spaces. The soft white finish adds a clean, professional look that complements any clinical interior.

##### **Dimensions:**

Trolley 60x41x89cm

Top shelf 53x29.5cm – middle 27.5x47.5cm – bottom 54x34cm