



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/dental-cart-soft-trolley-white-3-layers-600x410x890-mm-lxpxh/>



## Equipment trolley, white, 3 shelves, 600x410x890 mm (lxpxh)

The ES1003BS Dental Cart is a fundamental element for any medical practice or laboratory thanks to its versatility and user-friendliness. 600X410X890 MM, white, 3 layers

**361,32 €**

**Available**

**SKU:** ES1003BS

**Categories:** Multifunctional cart

### Description

---

Dental Cart SOFT 3 layers 600X410X890 MM, white, 3 layers

The ES1003BS Dental Cart is a fundamental element for any medical or dental practice, laboratory or beauty center, thanks to its versatility and user-friendliness.

#### Main features Dental Cart SOFT 3 layers

- Made of metal, highly resistant;
- easy to clean and sanitize;
- Smooth movement;
- 3 layers (up 53x29,5 cm, middle 27,5x47,5 cm, down 54x34 cm);
- middle layer adjustable in height;
- dimensions 60x41x89cm.