



Dentalcase mobile veterinary unit, 4 instruments and internal compressor

Dentalcase mobile veterinary unit with DCI style syringe and 2 cords, internal compressor

Available

SKU: UV4VA

Categories: Veterinary Dental Units

Description

Dentalcase mobile veterinary unit.

Minimum size, maximum efficiency! These mobile veterinary dental units permit unobstructed operations.

The power supply of the spray water tank, the vacuum suction system and the silent-running, built-in compressor make them fully independent and let you transport and position them anywhere you want.

Functions are controlled by simple and intuitive hand and foot controls. These units offer all the functions you need to provide effective veterinary dental treatments. They are robustly made from top quality materials and feature intelligent and stylish designs.

The front panel houses the water/ spray regulator and the instrument selectors. The bacteriological suction filter is positioned on the back of the trolley.

With its rolling luggage design, Dentacase makes you feel as if you are going on holiday!

Technical features Dentalcase mobile veterinary unit.

Power Supply	220/240 V – 50/60 Hz
Tank capacity	6 litres
Dimensions (w x h x l)	600 x 410 x 290 mm
Weight	30 Kg



DCI style syringe	<ul style="list-style-type: none">•
High speed handpiece /pneumatic micromotor cord	<ul style="list-style-type: none">• 2
Foot control	Disc type, progressive
Noise level	50 dB(A)
1lt spray water tank	<ul style="list-style-type: none">•
Suction system	Vacuum
Bacteriological suction filter	<ul style="list-style-type: none">•
Suction system	Vacuum
Warranty	24 months