



## Dentalcase mobile veterinary unit, 4 instruments and internal compressor

Dentalcase mobile veterinary unit with DCI style syringe and 2 cords, internal compressor

**Available**

**SKU:** UV4VA

**Categories:** Veterinary Dental Units

### Description

Dentalcase mobile veterinary unit.

Minimum size, maximum efficiency! These mobile veterinary dental units permit unobstructed operations.

The power supply of the spray water tank, the vacuum suction system and the silent-running, built-in compressor make them fully independent and let you transport and position them anywhere you want.

Functions are controlled by simple and intuitive hand and foot controls. These units offer all the functions you need to provide effective veterinary dental treatments. They are robustly made from top quality materials and feature intelligent and stylish designs.

The front panel houses the water/ spray regulator and the instrument selectors. The bacteriological suction filter is positioned on the back of the trolley.

With its rolling luggage design, Dentacase makes you feel as if you are going on holiday!

Technical features Dentalcase mobile veterinary unit.

|                        |                      |
|------------------------|----------------------|
| Power Supply           | 220/240 V – 50/60 Hz |
| Tank capacity          | 6 litres             |
| Dimensions (w x h x l) | 600 x 410 x 290 mm   |
| Weight                 | 30 Kg                |



|   |                        |
|---|------------------------|
| DCI style syringe                               | •                      |
| High speed handpiece /pneumatic micromotor cord | • 2                    |
| Foot control                                    | Disc type, progressive |
| Noise level                                     | 50 dB(A)               |
| 1lt spray water tank                            | •                      |
| Suction system                                  | Vacuum                 |
| Bacteriological suction filter                  | •                      |
| Suction system                                  | Vacuum                 |
| Warranty  | 24 months              |