



Dentalcase mobile veterinary unit, 4 instruments and internal compressor

Dentalcase mobile veterinary unit with DCI style syringe and 2 cords, internal compressor

Available

SKU: UV4VA

Categories: Veterinary Dental Units

Description

The **Dentalcase mobile veterinary dental unit** is designed to provide an efficient, autonomous and easy-to-transport workstation, ideal for managing veterinary dentistry procedures with minimum space requirements and maximum operating practicality.

Its **trolley structure** allows the unit to be moved wherever needed, supporting organised workflows without obstructions. This compact and functional solution is designed for professionals who need mobility without giving up the essential functions required during procedures.

The **spray water tank supply**, the **air vacuum suction system** and the **built-in silent compressor** make Dentalcase fully autonomous, allowing practical use in different operating environments.

The controls are **simple and intuitive**, and can be operated by both hand and foot. Front control panel houses the **water-spray regulators** and the **instrument lever selectors**, while the bacteriological suction filter is positioned on the rear of the trolley.

Made from quality materials, Dentalcase ensures strength, reliability and a refined design, providing all the necessary functions for veterinary procedures.

Technical features Dentalcase mobile veterinary unit.

Power Supply	220/240 V – 50/60 Hz
Tank capacity	6 litres
Dimensions (w x h x l)	600 x 410 x 290 mm
Weight	30 Kg



DCI style syringe	•
High speed handpiece /pneumatic micromotor cord	• 2
Foot control	Disc type, progressive
Noise level	50 dB(A)
1lt spray water tank	•
Suction system	Vacuum
Bacteriological suction filter	•
Suction system	Vacuum
Warranty	24 months