



Automatic sanitiser dispenser touch-free, 350 ml, 8 pz availables

Automatic sanitiser dispenser with infrared sensor for touch-free, fast and hygienic foam dispensing.

Available

SKU: DE2000B.7

Categories: Professional Equipment

Description

This automatic sanitiser dispenser is designed for quick, touch-free hand hygiene in professional environments. Equipped with an infrared sensor, it dispenses foam automatically when the hand is placed within the detection range, helping to reduce contact with surfaces and improve day-to-day hygiene management.

Its dispensing system converts liquid into foam, allowing smoother application and helping to reduce product waste by up to 50%. The anti-drip function helps keep the surrounding area cleaner, while the fast dispensing time of just 0.25 seconds ensures practical and efficient use.

Features automatic sanitiser dispenser

- Infrared sensor
- Touch-free automatic dispensing
- Liquid-to-foam conversion system
- Anti-drip function
- Waste reduction of up to 50%
- Detection range: 0-5cm
- Dispensing time: 0.25 seconds
- Dimensions: 73x73x190mm (LxDxH)
- Power supply: 4x1.5V AA batteries not included