



## Demineralizer, double

Double demineralizer DTPO2N system, compact size, designed to be directly connected to auto-load autoclaves.

**Available on request**

**SKU:** DTP02N

**Categories:** Demineralizers

### Description

The **Double Demineralizer System DTP02N** is an extremely compact universal device designed to remove all salts from water using special ion-exchange resins. It provides a continuous supply of pure water and is designed to be directly connected to automatic-load autoclaves.

When combined with the DS55 dispensing gun, it is the ideal optional accessory for rinsing instruments or filling manual-load autoclaves. The system is conveniently installed inside the sink cabinet, saving space and ensuring easy integration.

Tecnomed Italia universal demineralizers effectively remove dissolved salts from tap water. The ion-exchange resin columns retain the salts and deliver high-quality demineralized water.

The output water quality ranges from 0.1  $\mu$ S to 30  $\mu$ S, depending on the concentration of salts in the incoming water. This performance fully complies with EN13060 standards regarding water purity requirements for small sterilizers and beyond.

- 2 cartridge containers
- 2 empty cartridges
- 1 flow regulator
- 1 key
- 1 bracket
- 1 bucket containing resin (1,2 lt.)

Technical features double demineralizer DTPO2N system

Article Code	DTP02N
--------------	--------



Flow	48 l/h
Initial water purity	0.1 ?S
Max. water temperature	30°C
Internal cartridge size (d x h) mm	72,5 x 250
Total footprint (l x w x h) mm	270 x 145 x 335
Net weight	2,5 Kg
Assembly time	30 minutes
Power	230/12 vac
Warranty	24 months