



Surgical scrub sink stainless steel double station clinical lever tap

Wall-mounted **double station surgical scrub sink** in AISI 304 stainless steel, Scotch Brite finish, 1.5 mm thickness.

Designed for simultaneous use in operating theatres and high-traffic pre-operative scrub areas.

Equipped with dual clinical lever mixer taps and independent thermostatic mixing valves for precise temperature control.

Made in Italy.

Available on request

SKU: WASHCH-2

Categories: Sinks

Description

The wall-mounted double station surgical scrub sink is designed for clinical environments requiring simultaneous scrub access for two operators.

Manufactured in **AISI 304 stainless steel with Scotch Brite finish** and 1.5 mm thickness, it delivers structural robustness and long-term performance under intensive daily operation.

The extended basin geometry, with 10 mm radius rounded internal corners and a sloped base, ensures complete drainage and minimises water stagnation, supporting strict infection control protocols.

Each station operates independently and is equipped with a clinical lever mixer tap and a dedicated thermostatic mixing valve, allowing precise temperature regulation for enhanced safety and workflow efficiency.

The integrated high rear upstand and seamless front casing simplify sanitation procedures and reduce infiltration risks.

Designed and manufactured in Italy.

KEY FEATURES

- AISI 304 stainless steel structure, Scotch Brite finish (AISI 316 available upon request)
- 1.5 mm steel thickness
- Deep basin with 10 mm radius rounded internal corners



- Sloped base for complete drainage and reduced contamination risk
- Integrated high rear upstand to prevent water ingress
- Smooth front casing without joints or gaps
- Dual clinical lever mixer taps
- Dedicated thermostatic mixing valves for each station
- All technical components integrated within the upper panel, removable without tools
- Wall-mounted installation with dedicated suspension profiles (suitable for plasterboard walls)
- Stainless steel soap dispenser included
- Configuration: double station

DIMENSIONS

- Width (L): 1300 mm
- Depth (P): 590 mm
- Height (H): 1038 mm

Full dimensional drawings are available in the technical image provided.